Airtraq^{sp}

VIDEO LARYNGOSCOPE FOR DOUBLE LUMEN ENDOTRACHEAL TUBES

US Patent No 6,843,769
INSTRUCTIONS FOR USE
ENGLISH

DESCRIPTION AND INDICATIONS

The Airtraq SP DL is a SINGLE USE video laryngoscope designed to facilitate insertion of double lumen endobronchial tubes (DLT). It allows full visualization of the airway, during 100% of the intubation. It does not require hyper extension of the neck and permits intubating patients in virtually any position. It is provided clean, ready to use. Visualization can be performed directly through the eyecup or connecting it to any Endo Cam or to the accessories offered by the manufacturer.

SIZES

For use with double lumen tubes sizes 28-41 Fr.

Any style of DLT (right-sided or left-sided, Carina hooked or not hooked) can be used.

Minimum patient mouth opening: 18 mm.

TECHNIQUE FOR USING THE AIRTRAQ SP DL

I. PREPARATION

- Turn On the light. Switch located below the battery cover. If no light is seen, replace unit.

NOTE: The light stops blinking (30 sec) once the anti-fog system is fully activated.

- Remove the stylet from inside DLT.
- Check DLT cuffs and fully deflate both.
- Lubricate (preferably using water soluble lubricant, not spray) the lower 1/3 of the DLT and place it into the lateral channel of the Airtraq SP DL.

ORIENTATION OF THE DLT IN THE AIRTRAQ SP DL GUIDE CHANNEL

- Load the DLT and orient the deviated tip of the DLT so it points towards the light of the Airtraq SP DL. If the DLT tube is a right side stylet, then rotate the DLT 180 degrees until the deviated tip of the right DLT tube is pointing towards the light of the Airtraq SP DL.

- If the DLT has a carina hook, the hook should point toward the open side of the Airtraq SP DL guiding channel.

- Lubricate the DLT and place it into the lateral channel of the Airtraq SP DL without contacting the lens.

- Align the tip of the DLT with the end of the guide channel.

II. AIRTRAQ SP DL PLACEMENT (Fig. 1)

- Insert the Airtraq SP DL along midline of the patient's mouth. Take care to avoid pushing the tongue inside the oropharynx.
- Advance the Airtraq SP DL through the oropharynx keeping it in the midline.

- Before it reaches the vertical plane, begin looking to identify airway structures.
- Continue insertion until the epiglottis is identified. Place the tip of the Airtraq SP DL in the vallecula. Alternatively, the tip can be placed under the epiglottis, lifting it out of the way.

- Gently lift up the Airtraq SP DL to expose the vocal cords.

III. DLT INSERTION (Fig. 2)

- **Align the vocal cords in the center of the visual field** by gently moving the tip of the Airtraq SP DL as needed.

- Gently advance the DLT while in the guiding channel. STOP advancing when the DLT proximal cuff has just passed through the vocal cords. Do not insert further. This places the DLT at approximately the midpoint of the trachea. Check insertion depth.

IV. AIRTRAQ SP DL REMOVAL (Fig. 3)

- Separate the DLT from the Airtraq SP DL by pulling it laterally from the guiding channel, while holding the DLT in position.
- Check final insertion depth and location in the bronchi pursuant to hospital standard protocol, e.g. stethoscope or fiberoptic bronchoscope).

- Remove the Airtraq SP DL from the patient's airway following the midline.

WARNINGS AND PRECAUTIONS

- This product should only be used by personnel trained in insertion of Double Lumen Endobronchial tubes (DLT).

- Do not limitation.
- Do not incinerate unless batteries have been removed.
- Do not submerge in liquids.
- Use only with non-flammable anesthetics.
- Do not put pressure on the teeth with this device.
- Do not force the Airtraq SP DL into the upper airway.

USAGE TIPS

1. Initial experience should be gained in non-difficult airways.
2. Lubricate DLT generously (preferably using water soluble lubricants, not spray) and check that cuffs are fully deflated).
3. Insert the Airtraq SP DL, avoiding the tongue, and slide it softly and slowly.
4. Keep the Airtraq SP DL in the mouth's midline.
5. Look before the Airtraq SP DL gets to the vertical plane.
6. Do not insert too deep. If structures (arytenoids, epiglottis, etc.) are not clearly recognized or tube cannot be inserted, withdraw the Airtraq SP DL slightly.
7. Once the tip is located at the epiglottis, either at the vallecula (Macintosh style), or under the epiglottis (Miller style), gently lift up the Airtraq SP DL (do not tilt or use a lever action).
8. Advance the DLT slowly. If needed rotate DLT inside the channel.

BATTERY CHARACTERISTICS

Each Airtraq SP DL is equipped with two AAA batteries in serial connection that provide a voltage of 3 volts. The batteries provide power to the LED light and to the anti-fog system when the switch is turned on. They should not be replaced.

STORAGE, TRANSPORT, SHELF LIFE AND SERVICE LIFE

The Airtraq SP DL should not be used, stored or transported at temperatures below 0 °C / 32°F or over 40° C/104°F. The relative humidity must not exceed 95 %. The air pressure must not exceed 500 to 1060 hPa.

Airtraq SP DL shelf life is limited to the expiration date.

Airtraq SP DL service life is limited to 40 accumulative minutes, 5 minutes before its end; the Airtraq SP DL will start blinking for 2 seconds every 20 seconds. Once maximum operating time is reached the Airtraq SP DL will blink continuously until battery is depleted.

DISPOSAL

To dispose the Airtraq SP DL once it has been used:

- Separate the eyecup from the main body by pulling it up.
- Remove the battery cover by pulling it away from the main body (pull away from the small notches).
- Remove the batteries from the Airtraq SP DL and place them in an appropriate battery recycling container (dispose of them according to established recycling policies). The batteries are classified as non hazardous waste and must comply with European Directive WEEE. However, the manufacturer recommends separating them from standard trash.
- Discard the Airtraq SP DL as any other potentially contaminated waste.


MANUFACTURER'S WARRANTY

The manufacturer warrants the Airtraq SP DL against faulty materials or manufacturing defects for only one use or until the expiration date, whichever comes first, provided that the Airtraq SP DL is used in accordance with the procedures set forth in these instructions. This Warranty is applicable only if the device is purchased from an authorized distributor.

The Airtraq SP DL is designed for SINGLE-PATIENT-USE.

This device has not been designed to be cleaned or sterilized. Use beyond this recommendation may generate serious consequences in the product's performance and will void the Airtraq SP's warranty. The manufacturer disclaims all other warranties, whether expressed or implied, including,

without limitation, the warranties of merchantability or fitness for a particular use.

Manufactured by:

PRODOL MEDITEC LIMITED

No. 18, 7th Science Ave.,
Hi-Tec Coast, Zhuhai, Guangdong,
519085 P.R. China

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Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

For further advice on using the Airtraq SP please visit: www.airtraq.com or contact: info@airtraq.com or info.usa@airtraq.com
Airtraq is a registered trademark.

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