Airtraq^{sp}

VIDEO LARYNGOSCOPE

US Patent No 6,843,769
INSTRUCCIONES FOR USE

ENGLISH Orotracheal Intubation

DESCRIPTION AND INDICATIONS

The Airtraq SP is a SINGLE USE video laryngoscope designed to facilitate intubation. It allows full visualization of the airway, during 100% of the intubation. It does not require hyper extension of the neck and permits intubating patients in virtually any position. It is provided clean, ready to use. Visualization can be performed directly through the eyecup or connecting it to any Endo Cam or to the accessories offered by the manufacturer.

SIZES

Regular: A-011/ATQ-011

Blue. For use with ETT 7.0 – 8.5
Minimum patient mouth opening: 16 mm

Small: A-021/ATQ-021

Green. For use with ETT 6.0 – 7.5
Minimum patient mouth opening: 15 mm

Paediatric: A-031/ATQ-031

Purple. For use with ETT 4.0-5.5
Minimum patient mouth opening: 12 mm.

Infant: A-041/ATQ-041

Grey. For use with ETT 2.5-3.5
Minimum patient mouth opening: 11 mm.

TECHNIQUE FOR USING THE AIRTRAQ SP

I. PREPARATION

- Select the appropriate size Airtraq SP based on the size ETT to be used.
- Turn On the light. Switch located below the battery cover. If no light is seen, replace unit.

NOTE: The light stops blinking (30 sec) once the anti-fog system is fully activated.

- Lubricate the ETT and place it into the lateral channel of the Airtraq SP without contacting the lens.

- Align the tip of the ETT with the end of the lateral channel.

II. AIRTRAQ SP PLACEMENT (Fig. 1)

- Insert the Airtraq SP into the midline of the patient's mouth. Take special care to avoid pushing the tongue inside the oropharynx.

- Before it reaches the vertical plane, begin looking to identify airway structures.

- Continue insertion until the epiglottis is identified. Place the tip of the Airtraq SP in the vallecula. Alternatively, the tip can be placed under the epiglottis, lifting it out of the way.

- Gently lift up the Airtraq SP to expose the vocal cords.

III. ETT INSERTION (Fig. 2)

- Align the vocal cords in the center of the visual field** by gently moving the tip of the Airtraq SP as needed.

- Gently advance the ETT in the lateral channel. If needed rotate ETT inside the channel. Check insertion depth.

- Inflate the ETT cuff as normal and check for proper positioning.

IV. AIRTRAQ SP REMOVAL (Fig. 3)

- Separate the ETT from the Airtraq SP by pulling it laterally from the ETT, while holding the ETT in position.

- Remove the Airtraq SP from the patient's airway following the midline.

WARNINGS AND PRECAUTIONS

- This product should only be used by personnel trained in insertion of endotracheal tubes.
- Do not sterilize.
- Do not incinerate unless batteries have been removed.
- Do not submerge in liquids.
- Use only with non-flammable anesthetics.
- Do not put pressure on the teeth with this device.
- Do not force the Airtraq SP into the

upper airway.

USAGE TIPS

- Initial experience should be gained in non-difficult airways.
- Insert the Airtraq SP, avoiding the tongue, and slide it softly and slowly.
- Keep the Airtraq SP in the mouth's midline.
- Look before the Airtraq SP gets to the vertical plane.
- Do not insert too deep. If structures (arytenoids, epiglottis, etc.) are not clearly recognized, withdraw the Airtraq Avant slightly.
- Once the tip is located at the epiglottis, either at the vallecula (Macintosh style), or under the epiglottis (Miller style), gently lift up the Airtraq SP (do not tilt or use a lever action).
- Advance the ETT slowly. If needed rotate ETT inside the channel.

Each Airtraq SP is equipped with two AAA batteries in serial connection that provide a voltage of 3 volts. The batteries provide power to the LED light and to the anti-fog system when the switch is turned on. They should not be replaced.

BATTERY CHARACTERISTICS

The Airtraq SP should not be used, stored or transported at temperatures below 0 °C/32°F or over 40 °C/104°F. The relative humidity must not exceed 95 %. The air pressure must not exceed 500 to 1060 hPa.

STORAGE, TRANSPORT, SHELF LIFE AND SERVICE LIFE

The Airtraq SP shelf life is limited to the expiration date.

Airtraq SP service life is limited to 40 accumulative minutes, 5 minutes before its end; the Airtraq SP will start blinking for 2 seconds every 20 seconds. Once maximum operating time is reached the Airtraq SP will blink continuously until battery is depleted.

DISPOSAL

To dispose the Airtraq SP once it has been used:

- Separate the eyecup from the main body by pulling it up.
- Remove the battery cover by pulling it away from the main body (pull away from the small notches).
- Remove the batteries from the Airtraq SP and place them in an appropriate battery recycling container (dispose of them according to established recycling policies). The batteries are classified as non hazardous waste material and comply with European Directive WEEE. However, the manufacturer recommends separating them from standard trash.
- Discard the Airtraq SP as any other potentially contaminated waste.

MANUFACTURER'S WARRANTY

The manufacturer warrants the Airtraq SP against faulty materials or manufacturing defects for only one use or until the expiration date, whichever comes first, provided that the Airtraq SP is used in accordance with the procedures set forth in these instructions. This Warranty is applicable only if the device is purchased from an authorized distributor.

The Airtraq SP is designed for SINGLE-PATIENT-USE.

This device has not been designed to be cleaned or sterilized. Use beyond this recommendation may generate serious consequences in the product's performance and will void the Airtraq SP's warranty. The manufacturer disclaims all other warranties, whether expressed or implied, including, without limitation, the warranties of merchantability or fitness for a particular use.

Manufactured by:
PRODOL MEDITEC LIMITED
No. 18, 7th Science Ave.,
Hi-Tec Coast, Zhuhai, Guangdong,
519085 P.R. China

USA Representative:
AIRTRAQ LLC
24301 Woodstage Dr.
Bonita Springs, FL 34134-2958, USA

EU Representative:
PRODOL MEDITEC, S.A.
Muelle Tomás Olavari 5, 3º
48930 Las Arenas (Vizcaya) SPAIN
519085 P.R. China

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

For further advice on using the Airtraq SP please visit:
www.airtraq.com or contact:
info@airtraq.com or info.usa@airtraq.com

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Airtraq is a registered trademark.